

BaD EIPh beamline

Changing Gratings/Prefocusing

October 4, 2010

- Close the valve of the shutter
- Close the valve of the main experimental chamber

“Grating”

- In the program “Main ABS 2010 v1.vi” choose one of the following “Grating”:
 - AlMgF2** for energy range **4.6 - 13 eV** (requires **Si prefocusing**)
 - SiC** for energy range **14 - 18 eV** (requires **Pt prefocusing**)
 - Pt** for energy range **14 - 40 eV** (requires **Pt prefocusing**)
- Go to the monochromator chamber and insert the required grating **manually**, putting the marker of the feedthrough at the following position:
 - For **AlMgF2**: **161.5mm (70)**
 - For **SiC**: **80.06mm (40)**
 - For **Pt**: **-1.3mm (30)**

“Prefocusing Mirror”

- Depending on the grating chosen, decide on Prefocusing mirror considering that:
 - For **AlMgF2** grating you need **Silicon pref. mirror**
 - For **SiC** grating you need **Platinum pref. mirror**
 - For **Pt** grating you need **Platinum pref. mirror**
- If you need to change the Prefocusing Mirror do the following operations:
 - In the program “Main ABS 2010 v1.vi” click on the “Select” button of the “Prefocusing Operation” and choose the correct mirror.
 - Press the “Change Prefoc” button and wait till the pop-up disappears.
 - The “Actual Prefocusing Mirror” will be upgraded to the selected prefocusing mirror.
 - Manually check if the correct mirror has been inserted.

“Energy”

- e. In the “Energy” field write the value of photon energy (**Note: here the decimal separator is the comma**) you desire.

Option 1:

- a. Choose the “**Undulator Gap**” in “**Auto Movement**” for maximizing the photon flux automatically using the 1st harmonic (horizontal polarization).
- b. Press the “Set value - OK” button for changing both the photon energy and the gap of the undulator.

Option 2:

- a. Choose the “**Undulator Gap**” in “**Fixed**”
- b. Press the “Set value - OK” button for changing just the photon energy (the undulator gap will not be changed).
- c. Put the required undulator gap value (**Note: only here the decimal separator is the dot**) in “Gap Value” of the “Undulator Gap” panel and press “Set gap”.

Prefocusing re-alignment

- f. Go to the photodiode LTM and insert the photodiode **manually** (put the red markers equal, at about 20mm).
- g. Open the valve of the shutter.
- h. Press the “In line DAQ” button for acquiring the photodiode current.
- i. In the “Manual Movements” panel select the “Motor” “Prefoc. Pitch (M5)”, “Direction” “Forward” or “Back”, “number of pulses” of “1”
- j. Press the “Set” button several times to maximize the “Photodiode Current” shown in the graph.
- k. In the “Manual Movements” panel select the “Motor” “Prefoc. Roll (M4)”, “Direction” “Forward” or “Back”, “number of pulses” of “2”
- l. Press the “Set” button several times to maximize the “Photodiode Current” shown in the graph.
- m. Repeat the above steps **from i. to l.** to get the maximum photodiode current.
- n. Press again the “In line DAQ” to switch off the photodiode current acquisition.
- o. Extract manually the photodiode (red marker up to 45 mm)
- p. Open the valve of the main experimental chamber.